bewell[®]
connect

Shanghai

international holding
Corp. GmbH (Europe)
Eiffestraße 80, 20537
Hamburg Germany
Tel: 0049-40-2513175
Fax: 0049-40-255726

Made in China

Shenzhen Dongdixin
Technology Co., Ltd

No.3 Building Xilibaimang Xusheng
Industrial Estate, 518108 Nanshan,
Shenzhen China.

Tel.: 0086-755-27652316
E-mail: service@e-caretalk.com
www.e-caretalk.com

MANUFACTURED FOR

European market /
Marché Européen

VISIOMED GROUP SA

112, Avenue Kléber
75116 PARIS - France
contact@visiomed-lab.fr
Tel : +33 (8) 92 350 366
(0,40€/mn)
www.visiomed-lab.com
www.bewell-connect.fr

CE
0197

MyPeriTens

BW-PFX1

By Visiomed[®]

- A** Dispositif dans sa boîte de rangement - chargement / Device in its storage box - loading / Inrichting in zijn opbergdoos - het laden / Dispositivo en su caja de almacenamiento - carga / Vorrichtung in ihrem Aufbewahrungsbehälter - Aufladen
- B** Points de contact de chargement / Charging contact points / Laadcontactpunten / Puntos de contacto de carga / Ladekontaktstellen
- C1** Indicateur de lumière de charge / Charging light indicator / Opladerlampje / Indicador de luz de carga / Ladeanzeige
- C2** Prise de charge micro-USB / Micro-USB charging socket / Micro-USB-laadcontactdoos / Conector de carga Micro-USB / Micro-USB-Ladebuchse
- D1** Electrodes en acier inoxydable / Stainless steel electrode / Elektroden van roestvrij staal / Electrodo de acero inoxidable / Elektroden aus rostfreiem Stahl
- D2** Manchon en silicone / Silicone sleeve / Siliconenslang / Manga de silicona / Silikonhülle
- D3** Bouton d'alimentation / Power Button / Uit-knop / Botón de encendido / Einschaltknopf
- D3** / Led bleue: appareil en marche / Blue led: running device / Blaue led: het apparaat is ingeschakeld / Led azul: la unidad está encendida / Blaue led: Das Gerät ist auf
- D3** / Led verte: appareil chargé / Green led: loaded device / Groene led: inrichting loaded / Led verde: aparato cargado / Grüne led: gerät geladen
- C1** / Led orange: socle en charge / Orange led: charging stand / Oranje led: opladen stand / Led naranja: soporte de carga / Orange geführt: Ladeständer

beWellconnect

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MyPeriTens

CE
0197

Medical device class IIa, compliant with European Standards and with part 15 of the FCC rules.

TABLE OF CONTENTS

1. NOTES ON SAFETY	19
1.1. WARNING.....	19
1.2. CAUTION	20
1.3. GENERAL PRECAUTIONS	21
1.4. INSTALLATION AND USAGE PRECAUTIONS.....	22
2. INTENDED USE	22
3. FEATURES	22
4. USE.....	23
4.1. SUGGESTIONS BEFORE USE	23
4.1.1. INDICATOR DESCRIPTION	23
4.1.2. CHARGING	23
4.1.3. POWER ON/OFF.....	24
4.2. DOWNLOAD APPLICATION	24
5. HOW TO MAINTAIN AND STORE THE UNIT.....	24
5.1. CLEANING	24
5.2. MAINTENANCE.....	25
5.3. DISPOSAL	25
6. SPECIFICATIONS	25
7. TROUBLESHOOTING.....	27

Find out how to install and use the BewellConnect® app by visiting our website at www.bewell-connect.fr/application-bewellconnect

Hereby, Visiomed declares that the radio equipment type MyPeriTens is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address: <https://bewell-connect.com/en-gb/products-conformity/>

The manufacturer keeps the right to modify the specification and technical characteristics of the product without prior notice.

Dear Customer,

You have just purchased MyPeriTens from BewellConnect®. In order to use it optimally and efficiently, we recommend that you read these operating instructions carefully.

1. NOTES ON SAFETY

It is important that you read all the warnings and precautions included in this manual. They are intended to keep you safe, prevent injury and avoid a situation that could result in damage to the device.

1.1. WARNING

DO NOT USE THIS DEVICE WITH THE FOLLOWING MEDICAL DEVICES:

- Using this stimulator with other electronic medical devices may cause erroneous operation of those devices.
- Internally transplanted electronic medical devices, such as a pacemaker.
- Electronic life support equipment, such as respirators.
- Electronic medical devices attached to the body, such as electrocardiographs.

DO NOT USE THIS DEVICE UNDER THESE CONDITIONS:

- The device may cause lethal rhythm disturbances in certain susceptible individuals, consult your physician before using this device.
- For hygiene reasons, do not share this device with others. This device is a single user device: meaning only one person should use this device. There is a risk of infection, even if you clean the device.
- You are pregnant. The safety of electrical stimulation during pregnancy has not been established.
«British Journal of Anaesthesia»: Transcutaneous electrical nerve stimulation. Iain Jones MB ChB FRCA Mark I Johnson PhD. doi:10.1093/bjaceaccp/mkp021. 26 June, 2008
- You must be 6 weeks postpelvic surgery or vaginal childbirth to use this device.
- People must not be treated with this device.
- If you are less than 6 weeks postpelvic surgery or vaginal childbirth.
- If you have active pelvic cancer.
- If you have complete peripheral denervation of the pelvic floor, nerve damage or other diseases that reduce pain sensation in your genital area.
- If pain killers that reduce the feeling of pain in the vaginal area have been administered.
- During menstruation, abnormal vaginal discharge, vaginal or urinary tract inflammation or infection.
- While using a sanitary tampon or diaphragm.
- While the user is connected to high-frequency surgical equipment; it may cause burn injuries on the skin under the electrodes, as well as problems with the stimulator.

- Near of shortwave or microwave therapy equipment, since this may affect the output power of the stimulator.

- In case of fever.

People must not be treated with this device in case of:

- Extra-urethral incontinence (fistula, ectopic ureter),
- Overflow incontinence due to outflow obstacle,
- Serious retention of urine in the upper urinary tract,
- Symptoms of bladder infection, urinary or vaginal infections, localized lesions or other undiagnosed symptoms,
- Total/subtotal prolapsed uterus/vagina without doctor advice,

Paris. Postpartum pelvic floor muscle training and abdominal rehabilitation: Guidelines X. Deffieux, S. Vieillefoss and al. Doi: 10.1016/j.jgyn.2015.09.023

- Epilepsy,
- Irregular menstrual bleeding cycles.

DO NOT USE ON THESE INDIVIDUALS

- Children or infants, because the device has not been evaluated for pediatric use.

- This device is not designed for use by people (including children) with reduced physical, emotional or mental abilities or those lacking experience or knowledge unless used in the presence of someone responsible for their security, under supervision or following instructions about how to use the device.

DO NOT USE THIS DEVICE DURING THESE ACTIVITIES

- When in the bath or shower; penetrating water can impair function.

- While sleeping.

- When driving, operating machinery, or similar actions needing fine motor control; electrical stimulation may cause brief involuntary muscle movements.

1.2. CAUTION

- Do not use MyPeriTens probe for more than 30 minutes per session without medical advice.

- If tissue irritation should occur, treatment should be temporarily discontinued. If problems continue, contact your health care provider. Hypersensitivity can occur in isolated cases.

- Always turn the power off before removing or changing the location.

- If the stimulator is not functioning properly or you feel discomfort, immediately stop using the device.

- If you have suspected or diagnosed heart disease, you should follow precautions recommended by your physician.

- If you have suspected or diagnosed epilepsy, you should follow precautions recommended by your physician.

- Consult with your physician prior to using the device after a recent surgical procedure, as stimulation may disrupt the healing process.

- Use caution if you have a tendency to bleed internally, such as following an injury or fracture.

- To minimize the risk of infection, wash your hands thoroughly – preferably with anti-bacterial soap – before touching the probe.

- You should stop using the device and consult with your physician if you experience adverse reactions from the device.

1.3. GENERAL PRECAUTIONS

- This product should not encourage self-medication or adaptation of your treatment.

- Use this device only for the use for which it is designed as described in this manual. Do not use for any other purpose except for what it is intended for.

- Inspect the equipment prior to use, if you detect any problem on the device, do not use it and contact customer service.

- Do not use the device if it is damaged. The continuous use of a damaged unit may cause injury, or serious danger.

- The long-term effects of electrical stimulation are unknown.

- Follow the maintenance instructions specified in this manual.

- Before and after use, the part of the electrode in contact with the body and the silicone sleeve should be cleaned – please consult chapter 5.

- Keep unit away from children. The probe is a choking hazard and the USB cable could cause strangulation. The unit contains small pieces that may be swallowed. If swallowed immediately contact your physician.

- Use this device only with the silicone sleeve, adaptor and storage case recommended by the manufacturer.

- Do not freeze or microwave the device.

- Do not maintain or repair the device while the device is in use.

- Do not modify this device without authorization of the manufacturer.

- Use this device in an ambient temperature range of 5°C to 40°C / 41°F to 104°F and humidity range : 15% to 90%.

- Do not expose this device to extreme temperature conditions > 50°C to < -10°C / > 122°F to < 14°F.

- This device must always be placed in a clean and dry place in its storage box. Do not store this device to extreme humidity (below 0% RH or over 93% RH). Failing to do so may affect the performance of the stimulator.

- Do not expose this device to sunlight or water (unless for light cleaning, see section 5). Never immerse the device and the storage case in a liquid.

- Do not use this device outdoors.

- Never drop the device.

- Do not attempt to open the device. In case of problems, contact customer service.
- Do not disassemble the sensor or charging case.
- Do not expose this device to electric shocks.
- If you have any problems with this device, such as setting up, maintaining or using, please contact customer service.

1.4. INSTALLATION AND USAGE PRECAUTIONS

- Please move the transparent label from the metal part before using the device.
- Before use, check and ensure the main device and all accessories are intact and not damaged in any way, as use of the probe in this condition may increase likelihood of infection.
- If you feel any discomfort, stop treatment immediately.
- Place the storage case away from splashing water.
- Connect the main device and app correctly and securely.
- Always adjust the output intensity in the comfort level. If you feel uncomfortable, adjust the output intensity or stop treatment.
- To clean the device, please consult the chapter 5.
- Do not use cream, lubricating oil or other product which can damage the device. If vaginal environment is dry, apply a small amount of water-based lubricant to the inserting end of the probe.
- MyPeriTens can be reused, but should be cleaned before and after each use, and used by one individual only.
- After use, always put the device in the storage box. Do not store other objects in this box.

2. INTENDED USE

MyPeriTens is a wireless and connected pelvic floor electrostimulation device designed to treat stress, urge and/or mixed urinary incontinence in women. It allows to rehabilitate weakened pelvic floor muscles and to relieve pain of these pelvic floor muscles.

It may also be used by a health care provider to treat specific pelvic floor conditions by electrical stimulation.

3. FEATURES

This device uses electrical stimulation therapy, mainly for neuromuscular stimulation. Features are shown below:

- MyPeriTens device contains three functions: electrical stimulation (EMS and TENS), vibration massage, and vibration guided exercise (more info in App User Manual available online)
- EMS: Electrical Muscular Stimulation or NMES: NeuroMuscular Electrical Stimulation
- TENS: Transcutaneous Electrical Nerve Stimulation
- The Bluetooth is used to achieve wireless control and make it convenient to use with the mobile application (available on iOS and Android).

- The BewellConnect® app named MyPeriTens will work as a remote of the MyPeriTens device and also provide Kegel exercise programs (more information in the user manual of the application on the website www.bewell-connect.fr). The ergonomic design effectively prevents removal of the vaginal probe and rotation, to ensure that treatment is effective.
- The edge is wrapped with silicone to ensure a comfortable touch when treating.
- The waterproof level rating of MyPeriTens reaches IP56, to ensure that it is convenient to clean.

4. USE

4.1. SUGGESTIONS BEFORE USE

4.1.1. INDICATOR DESCRIPTION

MyPeriTens has two different indicator lights;

- One is the orange/blue light on the vaginal probe;
- The other one is the orange/green light on the charging case.

1. Function of indicator light on probe

- When the device is turned on without a Bluetooth connection, the blue light flashes
- When the device successfully connects with Bluetooth, the blue light stays on
- The orange light means low battery, it will flash for 10 seconds then device will turn off

2. Function of indicator light on the charging case

- Still under charging mode: Orange light
- Charging completed: Green light
- When probe is not positioned correctly, the green light will last for a long time and orange light will flash; you should replace the probe in the charging case and connect the USB cable again.

4.1.2. CHARGING

 DO NOT USE THE PROBE WHILE CHARGING.

Before first use, charge the device completely for approximately 3 hours, at room temperature. When the Lithium-ion battery is fully charged, it contains sufficient energy to use about four times for electrical stimulation treatment.

Notes:

- If the device is not fully charged when you start a treatment, the battery may run out during the treatment.
 - Treatment is not possible when the device is charging.
 - To ensure normal use, please fully charge before the first use.
 - If the indicator light turns orange, or the mobile application prompts that the battery is empty, please charge the device.
1. Place the main device in the storage case correctly.
 2. Insert the micro-USB plug of the adapter cable into the micro-USB socket of the storage case.

3. Insert the USB plug into the socket of the adapter and insert the adapter into a wall socket.
When charging, the indicator light of the storage case turns orange, and it will turn green when the battery is fully charged.
4. When charging is complete, remove the micro-USB plug from the micro-USB socket of the storage case.

4.1.3. POWER ON/OFF

1. Power on

Take the main device from the storage case and turn on the device by pressing the Power button. The blue light will flash, and the device will be in standby mode.

2. Power off

- Device will automatically power off if there is no Bluetooth connection within 5 minutes under standby mode or with no connection after 5 minutes without treatment.
- In the Power ON mode. Turn off the main device by holding the power button.
- When the orange light flashes to indicate low battery, the device will power off after 10 seconds.

Application will also display a pop-up to remind user of low battery during treatment session.

4.2. DOWNLOAD APPLICATION

The following information may be updated.

Download the MyPeriTens application to your smartphone or tablet via the App Store or Google Play. Make sure that the operating system is up to date (iOS 8.0 or later Or Android 4.3 or later).

We invite you to read the user manual of the application carefully for more information. (on the website www.bewell-connect.fr) about exercises and function of the device.

5. HOW TO MAINTAIN AND STORE THE UNIT

5.1. CLEANING

After finishing treatment, close the application, shutdown the device, then remove it. Before and after each use, clean the probe and silicone sleeve with clean water or mild soap water. Wipe them with dry and soft cloth.
Store the device in the storage case, always placed in a clean and dry place.

Never immerse the device and the storage case in a liquid.
Do not clean the device in an autoclave or dishwasher.
Do not clean with boiling water, corrosive cleaning agent, chemicals. Do not dry with abrasive materials (Risk of contamination and subsequent infection).
Do not share this device with another user.

5.2. MAINTENANCE

Follow the maintenance instructions specified in this manual.

1. Maintenance and all repairs should only be carried out by an authorized agency. The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorized people.
2. The user must not attempt any repairs to the device or any of its accessories. Please contact customer service for repair.
3. Opening of the equipment by unauthorized agencies is not allowed and will terminate any claim to warranty.
4. Check the device before each use for signs of wear and/or damage. Replace worn items as required.
5. To ensure the safety and effectiveness of the use, please use the main device and accessories provided by the manufacturer.
6. If you do not use the machine for a long time, please charge it every 2-3 months to ensure lithium battery performance.

5.3. DISPOSAL

Please dispose of the device in accordance with the directive 2012/19/EU –WEEE (Waste of Electronic Equipment). If you have any questions, please refer to the local authorities responsible of waste disposal. Note that the lithium battery in the probe poses a fire hazard and may explode if incinerated.

6. SPECIFICATIONS

Product name	BewellConnect® MyPeriTens
Model	BW-PFX1
Intended use	This device is used for the treatment of stress, urge and mixed urinary incontinence. It allows to rehabilitate weakened pelvic floor muscles and to relieve pain of these pelvic floor muscles.
Pelvic floor electrical stimulation output parameters	
Frequency range	2~100Hz
Pulse width range	50~350 us
Maximum output voltage	80 V at 500 ohm load
Output intensity level	0~80 levels
Treatment time	1~99 min
Pelvic floor muscle training exercises output parameters	
Contract time	1-60 s

Relax time	1-60 s
Cycle period	1-99
Preset programs	More than 16 programs
Vibrations massage output parameters	
Output intensity level	1-10 levels
Preset programs	More than 16 programs
Working noise	< 60dB (Measuring distance is 50 cm, back ground noise < 40db)
Classification	Classe IIa Applied part type BF IP56
Dimensions	Ø31 x 128.5 mm / Ø1.2 x 5 in
Weight	65 g / 0.14 lb
Rated input - adapter	100-240 V
Rated output - adapter	5 V-300 mA
Frequency - adapter	50-60 Hz
Wireless maximum separation distance	10 m (in an open environment)
Wireless recommended separation distance	3 m
Normal working conditions	Temperature: 5°C - 40°C / 41°F - 104°F Humidity: 0-93% RH Atmospheric pressure: 70-106Kpa
Storage and transportation conditions	Temperature: -10°C - 50°C / 14°F - 122°F Humidity: 0%-9% RH Atmospheric pressure: 70-106Kpa
Technology	Bluetooth 4.0 Smart Ready
Battery	3.7 V / 250mAh Li-ion
Frequency	2400-2483.5 MHz
Channel bandwidth	2 MHz (-20 dB)

Types of spread spectrum	FHSS
Modulation type	GFSK
Number of channels	40
Channel spacing	2 MHz
Channel frequency	0-39 Channel 2.402-2.480 GHz
Dwell time (if FHSS)	400 ms
Hopping rate (if FHSS)	1600 Hz
Maximum output powers	4 dBm
Accessories provided	Main device, charging box, user manual, USB cable, adapter provided

7. TROUBLESHOOTING

ERROR	POSSIBLE CAUSES	SOLUTION
LED light on the main device flashes orange, and then turns off automatically.	The battery is dead	Charge the device.
LED light on the main device flashing orange in working mode.	Low power	Battery power is at less than 25%. After use, charge the device.
Stimulation is uncomfortable.	Is the intensity too high or too low?	Adjust the output intensity.
	The electrode poor contact with vaginal of sujet	Use a lubricant.
	The electrode has poor contact with the vagina	Adjust the placement according to this user manual.

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EMC-DATA

EMC data

1. This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
2. Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
3. Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
4. Caution: this machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Refer to further guidance below regarding the EMC environment in which the device should be used.

GUIDANCE AND MANUFACTURE'S DECLARATION – ELECTROMAGNETIC EMISSION		
The BW-PFX1 is intended for use in the electromagnetic environment specified below. The customer of the user of the BW-PFX1 should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The BW-PFX1 use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The BW-PFX1 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

GUIDANCE AND MANUFACTURE'S DECLARATION – ELECTROMAGNETIC IMMUNITY			
The BW-PFX1 is intended for use in the electromagnetic environment specified below. The customer or the user of BW-PFX1 should assure that it is used in such an environment.			
Immunity test	EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2kV for power supply lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BW-PFX1 requires continued operation during power mains interruptions, it is recommended that the BW-PFX1 be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V RMS outside the ISM band, 6 V RMS in the ISM and amateur radio bands	Portable and mobile RF communications equipment should be used no closer to any part of the BW-PFX1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	10 V/m	
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BW-PFX1 is used exceeds the applicable RF compliance level above, the DEVICE should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the DEVICE. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE DEVICE .			
The BW-PFX1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BW-PFX1 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BW-PFX1 as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output of transmitter (w)	Separation distance according to frequency of transmitter (m)		
	150 KHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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EN: Visiomed* will repair or replace this product free of charge in the case of defective parts or manufacturing defects, in accordance with the conditions mentioned below as follows:

DURATION: 24 MONTHS RETURN TO WORKSHOP

LIMITS AND EXCLUSIONS: This guarantee concerns only the original final purchaser. A purchase invoice, or another proof of purchase, with this guarantee card will be required to obtain an after-sales service, in accordance with this guarantee. This guarantee card will not be extended to another person only the original final purchaser. This guarantee becomes void if the serial numbers on the product are modified, replaced, illegible, absent, or if repair has been carried out by a service not approved, including the user.

This guarantee covers only the defects of the material or parts, occurring during normal use of the product. It does not cover the damage caused during the transport of the apparatus, causes due to repairs being carried out by the distributor, by any modifications undertaken, any connection of equipment not approved by Visiomed*, or causes contrary to those written in the user manual or notice. Moreover, the present guarantee does not cover damage due to falls, bad handling, bad installations, damage by fire, floods, lightning, or any other natural disaster. This guarantee does not cover the packing of the material, the accessories, the defects caused by commercial exposure of the product, show room, sale space, demonstration etc... Normal maintenance, cleaning and the replacement of parts where wear is normal, are not covered by the terms of this guarantee. Visiomed* and its representatives and agents will not in any case be held responsible for any damage and consecutive damages due to the mishandling of this product. This guarantee is the only valid one at Visiomed*, any other guarantee (commercial guarantee) except this one will not be taken into account.

IMPORTANT: During the guarantee period if you are dissatisfied with the repairs of this product, please contact the Visiomed* customer service.

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EN: Disposal of Waste Equipment by Users in Private Households in the European Union.

For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.

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EN: BewellConnect® is a brand registered by Visiomed Group SA. The constructor keeps the right to modify the photographs and characteristics without previous notice. Registered design. The logos and trademarks are the property of their respective owners.

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© 07/2017 BW-PFX1