bewell
connect

MANUFACTURED FOR

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
BW-TSX - BW-TS1

By **Visiomed**

Stimulator BW-TS1
Rating: 5V $\approx$ 300mA
Output Max: 60mA at 500-1000 Ohm
Manufactured for: BEWELL CONNECT CORP
575 Boylston Street, Ste. 3W
Boston, MA 02116 - USA

- 1** On / off button / Boton de encendido / apagado / Bouton Marche / Arrêt
- 2** Status indicator / Indicador de estado / Indicateur de statut
- 3** Removable connector / Conector extraíble / Connecteur amovible
- 4** Micro-USB socket of the device / Conector micro-USB del dispositivo / Prise micro-USB de l'appareil
- 5** Electrodes pads / Electrodes almohadillas / Plaquettes d'électrodes
- 6** Micro-USB cable / Cable micro-USB / Câble micro-USB

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MyTens - MyTens Pro

This MyTens device complies with all the applicable standards and regulations related to exposure to electromagnetic fields.

The device complies with all the standards relating to Class II electrical medical devices and to devices which use electrical stimulation for use at home.

Electromagnetic compatibility (EMC) : The MyTens device complies with safety standards EMC IEC 60601-1-2. It is designed to be used at home.

CONTENTS

1. INDICATIONS FOR USE	7
2. IMPORTANT SAFETY PRECAUTIONS AND WARNINGS	7
2.1 CONTRAINDICATION	7
2.2 WARNING	7
2.3 CAUTION	9
2.4 GENERAL PRECAUTIONS	10
3. CHARACTERISTICS	11
4. DIRECTIONS BEFORE USE	11
4.1 CHARGING	11
4.2 ELECTRODES	12
5. USE	13
5.1 PROGRAMS PRESENTATION	13
5.1.1. FORMS OF WAVES	18
5.1.2. TENS PROGRAMS	18
5.1.3. ELECTRICAL MUSCLE STIMULATION (EMS) PROGRAMS	20
5.2. DOWNLOADING THE MOBILE APPLICATION	23
5.3. STARTING AND USING MYTENS / MYTENS PRO	23
6. ACCESSORIES	23
7. STORAGE AND MAINTENANCE	24
7.1 CLEANING	24
9. TECHNICAL CHARACTERISTICS	25
10. DESCRIPTION OF THE WIRELESS FUNCTION	26
11. IMPORTANT INFORMATION REGARDING ELECTRO MAGNETIC COMPATIBILITY (EMC)	27

Find out how to install and use the BewellConnect® app by visiting our website at www.bewell-connect.us/bewellconnect-app.

The manufacturer reserves the right to change the product's technical specifications without prior notice.

Dear Customer,

You have just purchased Mytens from BewellConnect®. In order to use it optimally and efficiently, we recommend that you read these operating instructions carefully.

For specific information about vein, muscle, or other problems, PLEASE CONTACT YOUR DOCTOR.

The basic performance of the device is free from the production of unwanted or excessive stimulation output (device can only be used with our preset programs).

1. INDICATIONS FOR USE

TENS: The device is designed to be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, neck, upper extremities (arm) and lower extremities (leg) due to strain from exercise or normal household work activities.

EMS: The device is designed to be used for stimulate healthy muscles in order to improve and facilitate muscle performance.

2. IMPORTANT SAFETY PRECAUTIONS AND WARNINGS

2.1 CONTRAINDICATION

This stimulator must not be used in combination with the following medical devices:


- Internally transplanted electronic medical devices, such as a pacemaker.
- Electronic life support equipment, such as respirators.
- Electronic medical devices attached to the body, such as electrocardiographs.
- Using this stimulator with other electronic medical devices may cause erroneous operation of those devices.

2.2 WARNING

- Please read this manual carefully and always follow the treatment instructions.
- Only use this device for its intended purpose as described in these instructions.
- Consult with your physician before using this device, because the device may cause lethal rhythm disturbances to the heart in susceptible individuals.
- This device must always be placed in a clean and dry place.
- Do not expose this device to sunlight or water.
- Do not expose this device to electric shocks.
- Never drop the device.
- Follow the maintenance instructions specified in this manual.
- Do not attempt to open the device. In case of problems, contact customer service.
- This is a medical device. Keep out of the reach of children.
- Discontinue use of the device in case of anomalies or malfunction.
- Do not apply stimulation in the presence of electronic monitoring equipment (e.g., cardiac monitors, ECG alarms), which may not operate properly when the electrical stimulation device is in use.
- Stimulation should not take place while the user is connected to high-frequency surgical equipment, it may cause burn injuries on the skin under the electrodes, as well as problems with the stimulator.
- Do not use the stimulator in the vicinity of shortwave or microwave therapy equipment, since this may affect the output power of the stimulator.

- Possible effects from RF sources in the vicinity of the device (e.g, electromagnetic security systems, cellular telephones, RFID or other in- band transmitters).
- Do not operate this unit in an environment where other devices are used that intentionally radiates electromagnetic energy in an unshielded manner, such as MRI, RFID, metal detectors, and electronic article surveillance (EAS) anti-theft systems.
- Do not use the device if it is closer than 12 in (30.5 cm) to wireless communication equipment, such as wireless home network routers, cell phones, cordless phones and their base stations, and walkie-talkies. The electromagnetic interference of this wireless communication equipment may prevent the device from operating properly.
- This device is not designed to be used by persons (including children) whose physical, sensory or mental capabilities are reduced, or persons without experience or knowledge, unless they have been able to benefit, by the intermediary of a person responsible for their safety, from supervision or prior instructions concerning use of the device. It is possible that they will not be able to use it in accordance with the instructions of this user manual and be disturbed by the treatment.
- Do not modify the device or the electrodes. This could cause a malfunction.

Do not treat the the following parts or areas:

	Do not apply stimulation on your chest, the introduction of electric current on this area can cause heart rhythm disturbances, with a risk of death.
	Its effects on the brain are unknown. Stimulation to the head or on each side of the skull must be avoided.
	Stimulation on the sides of the neck or on the carotid artery can cause serious adverse effects on your blood pressure or your heart rhythm.
	Stimulation on the front of the neck can cause severe muscle spasms that can block your airways and cause breathing difficulties.
	Do not place the electrodes on the inside of body cavities, such as in the mouth. This device is only designed for external application.
	Do not place the electrodes on the genitals. They could stimulate inappropriate muscles or organs.

- Do not position the electrodes on broken or injured skin, or that which is dirty or unhealthy. Skin with irritation, sores or other lesions can lead to the injection of too much current on the area, which can cause burns.
- Do not place the electrodes near cancerous lesions because this may have a negative impact on these injuries.
- Do not place the electrodes on skin areas on which sensations are not normal. You may burn yourself due to a lack of perception of the high intensity of the current.
- Do not put the electrodes on areas that are swollen, red, infected or inflamed or on skin rashes (eg. phlebitis, thrombophlebitis and varicose veins). Stimulation should not be performed on

areas of thrombosis or thrombophlebitis because it can promote the circulation and lead to a greater risk of embolism.

- Do not put the electrodes on redness or open wounds. Open wounds may lead to applying too much current on the area, causing burns. They can also further the penetration of substances from the electrode into the skin.
 - Do not make sudden movements during a session. This could cause a dysfunction of the device.
 - If batteries leak and come into contact with the skin or eyes, wash immediately with copious amounts of water.
 - Do not apply stimulation over, or in proximity to, cancerous lesions.
 - Apply stimulation only to normal, intact, clean, healthy skin.
- Do not use the device in the following conditions:
- Do not use the device to children or infants, because the device has not been evaluated for pediatric use.
 - Do not use the device to pregnant women, because the safety of electrical stimulation during pregnancy has not been established.
 - Do not use the device to persons incapable of expressing their thoughts or intentions.
 - Do not use the electro-stimulator if you are monitored by a doctor and you have not consulted him before using it.
 - In the case of internal bleeding due to impacts or injury, do not use the device.
 - Do not use the device while driving, operating machinery or any other activity during which the electrical stimulation may lead to a risk of injury.
 - Do not use the device if you are subject to falling asleep during the session, as this may cause you to feel pain too late. If using at the time of going to bed, set the timer so that the device does not turn itself off automatically.
 - Never use MyTens in contact with water (in the bathroom, in the shower or in the pool, etc.) because this increases the risk of an electric shock and skin burns.

PAIN MANAGEMENT WARNINGS:

- If your pain does not improve, becomes seriously chronic or severe, or continues for more than five days, stop using the device and consult with your physician.
- The mere existence of pain functions as a very important warning telling us that something is wrong. Therefore, if you suffer from any serious illness, consult your physician in order to confirm that it is advisable for you to use this device.

WARNINGS AND PRECAUTIONS REGARDING THE ELECTRODE PADS:

- If you experience any skin irritation or redness after a session, do not continue stimulation in that area of the skin.
- Electrode pads should not touch each other when placed onto your skin.

2.3 CAUTION

- Do not bend or fold because the electrode pad may not function properly. Place the electrode pads onto the plastic film and then store into the sealed package when not in use.
- Do not apply ointment or any solvent to the electrode pads or to your skin because it will disrupt the electrode pads from functioning properly.

- Always place the electrodes in accordance with the electrode illustrations.
- Do not place on your spine or backbone.
- Electrode pad should not touch any metal object, such as a belt buckle or necklace.
- Do not share electrode pads with another person. This may cause a skin irritation or infection. Electrode pads are intended for use by one person.
- Always turn the power off before removing or changing the electrode pad location.
- If the stimulator is not functioning properly or you feel discomfort, immediately stop using the device.
- Dispose of the device, batteries, and components according to applicable legal regulations. Unlawful disposal may cause environmental pollution.

The effectiveness of your sessions will depend in part on how you place your electrodes on the skin. Before placing the electrodes, ensure you have good sensation in the area. If in doubt about the ideal placement you can move the electrodes to find the best sensation. If you find an area to be effective, note the placement for future use. If you suffer discomfort (unpleasant contractions during burst electrical neurostimulation programs), it is advisable to move the electrodes. If during electrical muscle stimulation programs (intended especially to create a contraction of the muscle) no contraction occurs, this may mean an intensity is too low or the electrodes have been placed incorrectly.

2.4 GENERAL PRECAUTIONS

- Before using the device, consult your doctor in the case: you take insulin for diabetes; muscle contractions stimulate glucose uptake / recent surgery; the use of this device may affect the healing process.
- It is necessary to complete the treatment in the MyTens application or press the on/off button of the device before removing the electrodes. If you do not complete the treatment and remove the electrodes or the device, and if you touch the connectors, it is possible to feel pain in the fingers. This sensation, however, is not dangerous.
- Before any use of electrodes, check their expiration date on the packaging. If the date has expired, do not use them.
- Use new electrodes when the electrodes are damaged, dirty, less adhesive or if you begin to feel discomfort during stimulation (uncomfortable pinching, tingling).
- This product must only be used with the adapter, the cable and accessories recommended by the manufacturer.
- Before use, check that the device and the electrodes are not damaged. In case of any damage, do not use the electrodes.
- Usable indoors and outdoors, the device is not designed to withstand all weather conditions.
- The device is not waterproof. Do not get wet and never use it in a damp environment.
- Do not use plaster or tape to attach the electrodes to the skin.
- For storing the electrodes, please follow the instructions in section 5.1 "Storage of the electrodes".
- In case of abnormal operation of the device in accordance with this user manual, discontinue use and contact customer service. Consult section 6 "Troubleshooting Guide" for information concerning potential malfunctions of the device.
- The patient is intend operator.

- The long-term effects of electrical stimulation are unknown.
- TENS is not effective in treating the original source or cause of the pain, including headache.
- TENS is not a substitute for pain medications and other pain management therapies.
- TENS is a symptomatic treatment and, as such, suppresses the sensation of pain that would otherwise serve as a protective mechanism.
- If you have suspected or diagnosed heart disease, you should follow precautions recommended by your physician.
- If you have suspected or diagnosed epilepsy, you should follow precautions recommended by your physician.
- Use this device only with the electrodes, adapter and Micro-USB cable recommended by the manufacturer.
- Use the transcutaneous electrical neurostimulation programs in bursts (TENS Burst) at rest. They can cause a muscle contraction in the treatment area.
- The safety of electrical stimulation during pregnancy has not been established.
- Consult with your physician prior to using the device after a recent surgical procedure, because stimulation may disrupt the healing process.

Side effects

Side effects are few and often represented by skin irritation (redness) associated to electrodes. They usually disappear within two hours following a session of electro -stimulation. If they persist for more than 24 hours, consult your doctor.

Users can also present myoclonus (involuntary muscle contractions) if the stimulation intensity is too high, or sensory disturbances.

In case of side effects during the use such as nausea or fainting, burns at the locations of the electrodes, painful sensations or headache, stop using the device and consult your doctor.

3. CHARACTERISTICS

- Pain relief and healing of muscle tension (TENS: 11 preset programs)
- Strengthening of the muscular mass (EMS: 8 standard programs)
- Designed to be used for temporary relief of pain.
- Designed to be used for the stimulation of healthy muscles.

4. DIRECTIONS BEFORE USE

4.1 CHARGING

Charge the device before the first use. In normal conditions (ambient temperature), it takes 5 hours to charge the device. A full charge allows a battery life of about 8 hours continuously when functioning in normal conditions. Please charge the device fully before beginning a session. If the charge is not complete at the beginning of a program, the battery may become depleted before the end of the session. You cannot use the device when it is charging.

The life of a rechargeable battery depends on the number of recharging/run-down cycles it undergoes and how these cycles are performed. The service life of rechargeable battery is more than 300 recharging/run-down cycles. We provide the following suggestions for the longer life of the battery:

- 1) Whenever the device is not used frequently, we recommend recharging the battery once a month.

- 2) The device can only be used safely with the original recharger, DO NOT reassemble or change the specification of the recharger. So please always use the recharger which supply by manufacturer or distributor.

Plugging in the device figure B on page 2

1. Disconnect the cable used to connect the two connector units from the unit with the battery light and micro-USB socket.
2. Insert the micro-USB cable into the micro-USB socket for the respective connector unit.
3. Insert the USB plug into a computer to charge the device.
4. When the battery light goes into continuous orange (which means that the battery is being charged) to a continuous green (the battery is full), remove the micro-USB plug from the micro-USB socket of the connector unit with the battery indicator (see the chapter 4). Use the USB plug only for the use described above (to charge and connect to the cable connecting the two connectors units).
5. Reconnect the cable connecting the two connectors units : Insert the unit with the cable into the micro-USB socket for the unit with the battery light .

4.2 ELECTRODES

BewellConnect® electrodes pads are used with MyTens devices, for muscle stimulation or pain relief. The electrodes pads help connect the MyTens muscle stimulator to your skin. Electrodes pads are the points where the pulses are transmitted to your motor nerves, and ultimately your muscle fibers.

- Ensure that the skin is healthy and shows no: irritated skin, injuries, open wounds, redness, rashes or pre-cancerous lesions.
- Position the electrodes only on intact and healthy muscles. In case of doubt, contact your doctor or physiotherapist.
- Clean skin thoroughly prior to each application of the electrodes. Electrodes do not stick well if any lotion, oil, make-up, dirt, etc. is left on skin. Always moisten skin with a damp cloth prior to applying the electrodes. This will improve the electrode adhesion and life.
- Remove the electrodes from the bag, ensure they are not damaged or worn. If the expiration date has past / is exceeded, do not use these electrodes (the life time for Mytens pads is two years).
- A BewellConnect® electrode pad is coated with a strong adhesive, developed to stick to your skin.
- Used in normal circumstances, a BewellConnect® electrode pad should last 15 to 20 sessions
- Remove the clear plastic film from the back of the electrode pad.
- Place electrodes on clean, dry and healthy skin near, and do not let them touch. Make sure there is a linear path between the two electrodes.
- Avoid any contact between both of the electrodes, do not place them on top of each other.
- Follow instruction ad steps within the app.
- Attach the electrodes to the device using the snap buttons provided for this purpose. (see

picture 5 page 2)

- Positioning the electrodes on the body is explained within the app for each program.
- It is possible that you may need additional assistance for placing the electrodes on certain areas (back area).
- Make sure that the device is switched off when moving the electrodes to avoid risk of injury to yourself or others.

Electrodes positions are displayed in the app. (Check the MyTens app user manual online)

WHEN SHOULD YOU REPLACE YOUR ELECTRODES PADS ?

- 1 When they no longer stick to your skin.
- 2 When the same intensity setting of the program creates less intense or uncomfortable muscle contractions.
- 3 When the adhesive is filled with lint, dust or hair.

5. USE

5.1 PROGRAMS PRESENTATION

INDEX PROGRAM ID	SUGGESTED TREATMENT SITE	TYPE OF ELEC-TRICAL STIMULATION	DESCRIPTIONS	CATEGORY OF PROGRAMS (IN THE APP)
1	SHOULDER	TENS ENDORPHINIC	This Pain relief strong massage program, use for temporary relief of pain associated with sore and aching muscles in shoulder, waist, back, neck, arms, or legs to strain from exercise or normal activities. Stimulate your circulatory system to help eliminate waste and toxins from your body. This should also induce an endorphin release , which react in the brain to reduce our perception of pain. This program is also good for back pain . Intensity should induce visible muscular contractions, but it must remain bearable so as not to cause any pain.	PAIN RELIEF + RELAXATION

2	LUMBAR	TENS MODULATED FREQUENCY	If you use the pain relief strong massage program often your body may become accustomed to the program, making it less effective. Use this Pain relief 1 program in combination or as an alternative for pain relief associated with sore and aching muscles in the shoulder, waist, back, neck, arms, or legs due to strain from exercise or normal activities. Intensity should be set to feel a tingling sensation; if it is too strong the program will become less efficient. You should not see your muscles moving.	PAIN RELIEF
3	ELBOW	TENS HANS	Use for management of pain to relieve symptoms from joint pain and arthritis. TENS HANS uses tiny electrical impulses sent through the skin to nerves to modify your pain perception. Intensity should be set to feel a tingling sensation; if it is too strong the program will become less efficient. You should not see your muscles moving.	PAIN RELIEF
4	FOREARM	TENS HANS	Use for management of pain to relieve symptoms from joint pain . TENS HANS uses tiny electrical impulses sent through the skin to nerves to modify your pain perception. Intensity should be set to feel a tingling sensation; if it is too strong the program will become less efficient. You should not see your muscles moving.	PAIN RELIEF
5	WRIST	TENS BURST	Use for temporary relief and management of pain associated with sore and aching muscles in the shoulder, waist, back, neck, arms, or legs due to strain from exercise or normal activities. Intensity should be set to feel a tingling sensation; if it is too strong the program will become less efficient. You should not see your muscles moving and feel muscular relaxation .	PAIN RELIEF + RELAXATION

6	HIP	TENS ENDORPHINIC	Use for muscle fatigue after a workout or training session. Stimulation will increase the blood flow to drain away toxins and improve and facilitate muscle performance. Intensity should be kept at comfortable level.	RELAXATION
7	THIGH	EMS 100 HZ	Allow muscles to exercise passively by sending electronic pulses to the muscle needing treatment. Improve muscle fatigue by recovering muscle strength and relieving tension. Intensity should be set at maximum bearable level, Intensity should induce light and bearable muscular contractions.	MUSCLE STIMULATION
8	KNEE	TENS GATE CONTROL MODULATED	If you use pain relief 1 program often your body may become accustomed to the program, making it less effective. Use this the pain relief level 2 in combination or as an alternative for pain relief associated with sore and aching muscles in the shoulder, waist, back, neck, arms, or legs due to strain from exercise or normal activities. Intensity should be set to feel a tingling sensation; if it is too strong the program will become less efficient.	PAIN RELIEF
9	LOW LEG	TENS GATE CONTROL MODULATED	If you use pain relief 1 program often your body may become accustomed to the program, making it less effective. Use this the pain relief level 2 in combination or as an alternative for pain relief associated with sore and aching muscles in the shoulder, waist, back, neck, arms, or legs due to strain from exercise or normal activities. Intensity should be set to feel a tingling sensation; if it is too strong the program will become less efficient.	PAIN RELIEF


10	ANKLE	TENS CON-VENTIONAL GATE CONTROL 100 HZ	Conventional Tens uses tiny electrical impulses sent through the skin to nerves to modify your pain perception. Use to help control and temporarily relieve specific or acute pain. The intensity must be increased gradually until there is tingling sensation that is pronounced without being painful. Electrodes should be positioned such a way to cover or surround the painful area	PAIN RELIEF
11	FOOT	TENS MODU-LATED FRE-QUENCY	Reduce pain and relieve muscle tension with this relaxing Comfort Massage program. Let your stress and tension melt away. Intensity should be gentle and tolerable	RELAXATION
20	SHOULDER	EMS MODU-LATED FRE-QUENCY	This Endurance program will induce contractions to help to condition the muscle in order to facilitate performance. Stimulation will increase the blood flow and the oxygen supply for faster recovery from muscle fatigue. The maximum tolerable stimulation intensity. The higher the stimulation intensity, the higher the number of muscle fibers being used.	MUSCLE STIMULATION
21	ABDOMEN	EMS MODU-LATED FRE-QUENCY	Use the Firming-up program regularly to help strengthen your muscles and recover muscle strength. The higher the intensity, the more effective the practice	FITNESS
22	LUMBAR	EMS MODU-LATED FRE-QUENCY	If your muscles are already firm, use this Reinforcement program to stimulate healthy muscles in order to improve performance. The maximum tolerable stimulation intensity. The higher the stimulation intensity, the higher the number of muscle fibers being used.	FITNESS + MUSCLE STIMULATION

23	FOREARM	EMS MODU-LATED FRE-QUENCY	Try this Workout prep program before a workout. Stimulation will increase the blood flow and the oxygen supply getting the muscle ready to work and contraction helps to condition the muscle to facilitate performance. Intensity should induce visible muscular contractions, but it must remain bearable. Use the maximum stimulation intensity. Allow one or two sessions to become accustomed to the program and gradually increasing the stimulation intensity every 3 or 4 contractions.	FITNESS
24	HIP	EMS MODU-LATED FRE-QUENCY	Stimulate healthy muscles to help recover muscle strength with this Strengthening level 1 program. Intensity should be set at maximum bearable level, particularly during muscular spasms.	FITNESS + MUSCLE STIMULATION
25	THIGH	EMS MODU-LATED FRE-QUENCY	Speed up muscle recovery; this Strengthening level 2 program displays contraction helps to condition the muscle to facilitate performance and regain muscle strength. Intensity should be set at maximum bearable level, particularly during muscular spasms.	FITNESS + MUSCLE STIMULATION
26	LOW LEG	EMS MODU-LATED FRE-QUENCY	Speed up muscle recovery; this Strengthening level 2 program displays contraction helps to condition the muscle to facilitate performance and regain muscle strength. Intensity should be set at maximum bearable level, particularly during muscular spasms.	FITNESS + MUSCLE STIMULATION

27	FOOT	EMS 50 HZ	Provides stimulation to facilitate recovery from muscle fatigue and to help recover muscle strength by sending electronic pulses to the muscle needing treatment; this causes the muscle to exercise passively. Intensity should induce light and bearable muscular contractions	FITNESS + MUSCLE STIMULATION
----	------	-----------	---	-------------------------------------


5.1.1. FORMS OF WAVES

*CONT: Form of continuous wave


*Cont = Continuous


BURST: Wave form deviation

HANS: Form of alternating wave

EMS: wave form

5.1.2. TENS PROGRAMS

Traditional transcutaneous neurostimulation program:

A common transcutaneous neurostimulation program is performed with a high frequency stimulation. This stimulation is intended to block the transmission of the pain signal to the nervous system. Generally, people use traditional transcutaneous neurostimulation for a session of about 30 minutes. It is possible to use these programs for a variable period and frequency according to what you want. It is normal to feel the effects of this program more intensely during the stimulation and to experience a decline after the end of the session. To improve effectiveness, adjust the intensity of the program, avoiding any painful or unpleasant sensation. Burst transcutaneous neurostimulation program:

A burst program (Burst TENS) allows low frequency stimulation in the form of a series of impulses. It enables the release of endorphins. It is possible to use the transcutaneous neurostimulation in bursts for 30 minutes. The frequency can be up to several times per day; for example, 3 times.

For a sensation of strong stimulation, it is advisable to adjust the intensity of the burst program. An unpleasant sensation may then appear. You can also adjust the intensity of the program at a level leading to perceptible muscle contractions in the stimulation area. These contractions are not dangerous. However, if you experience discomfort, it is advisable to decrease the intensity of the program, change or move the electrodes. The burst transcutaneous electrical neurostimulation (Burst TENS) can trigger a muscle contraction in the area of stimulation, so you should to perform this program at rest.

SUMMARY TABLE OF THE TENS PROGRAMS

TENS			IMPULSES			
Index program ID	Treatment time	Treatment phases	Frequency (time)	Width in µS	Wave form	Description
1	30min	Phase 1	2Hz (10s)	250	CONT	Continuous output
		Phase 2	4Hz (8s)	250		
		Phase 3	6Hz (6s)	250		
2	30min	Phase 1	80Hz (20s)	250	CONT	Continuous output
		Phase 2	80Hz (20s)	250		
		Phase 3	75Hz (4s)	250		
		Phase 4	10Hz (20s)	250		
		Phase 5	70Hz (4s)	250		
		Phase 6	65Hz (4s)	250		
3	30min	Phase 1	2Hz (3s)	150	HANS	Alternating output
		Phase 2	100Hz (3s)	200		
4	30min	Phase 1	2Hz (3s)	150	HANS	Alternating output
		Phase 2	100Hz (3s)	200	HANS	
5	30min	Phase 1	2Hz	250	CONT	Continuous output
6	30min	Phase 1	6Hz (30s)	250	CONT	Continuous output
		Phase 2	6Hz (30s)	250		
		Phase 3	8Hz (20s)	250		
		Phase 4	8Hz (20s)	250		
		Phase 5	10Hz (20s)	250		
		Phase 6	10Hz (20s)	250		
7	30min	Phase 1	100Hz	200	EMS	Up: 0.5s Cont: 7s Down: 0.5s Res : 7s

8	30min	Phase 1	100Hz	150	BURST	On: 0.25s Off: 0.25s
9	30min	Phase 1	100Hz	150	BURST	On: 0.25s Off: 0.25s
10	30min	Phase 1	100Hz	200	CONT	Continuous output
11	30min	Phase 1	40Hz (5s)	250	CONT	Continuous output
		Phase 2	6Hz (10s)	250		
		Phase 3	50Hz (5s)	250		

5.1.3. ELECTRICAL MUSCLE STIMULATION (EMS) PROGRAMS

EMS: Electrical Muscular Stimulation or NMES: NeuroMuscular Electrical Stimulation TENS: Transcutaneous Electrical Nerve Stimulation, emit impulses to cause a muscle contraction. They are especially used for strengthening the muscle mass. Only perform the electrical muscle stimulation program on intact and healthy muscles. Because they involve muscle contraction, they need to be performed at rest (lying or sitting).

When using for the first time, the electrical stimulation can lead to an atypical and strange sensation. It is advisable to begin the program on low intensity and change to a higher one to become used to it. The stimulation must remain comfortable and muscle contractions must not be harsh or painful.

Find more information about electrical muscle stimulation programs online (MyTens App User Manual and within the app MyTens available on Google Store and AppStore).

SUMMARY TABLE OF THE EMS PROGRAMS

EMS			IMPULSES			
Index program ID	Treatment time	Treatment phases	Frequency (time)	Width in μ S	Wave form	Description
20	28min	Phase 1	5Hz (5mn)	280	EMS	Preparation phase
		Phase 2	55Hz (10s)	280		Intensity: 50% Up: 1.5s Down: 0.75s
		Phase 3	6Hz (8s)	280		Intensity: 25%
		Phase 4	3Hz (10mn)	280		Phase of relaxation Intensity: 40%
21	22min	Phase 1	30Hz (2mn)	200	EMS	Up: 10s Cont: 5s Down: 10s Rest: 5s
		Phase 2	45Hz (9mn)	200		Up: 5s Cont: 5s Down: 5s Rest: 5s
		Phase 3	60Hz (9mn)	200		Up: 6s Cont: 8s Down: 6s Rest: 5s
		Phase 4	30Hz (2mn)	200		Up: 10s Cont: 5s Down: 10s Rest: 5s

22	28min	Phase 1	5Hz (5min)	300	EMS	Preparation phase
		Phase 2	55Hz (10s)	300		Intensity: 50% Up:1,5s Down:0,75s
		Phase 3	6Hz (8s)	300		Intensity: 25%
		Phase 4	3Hz (10min)	300		Phase of relaxation intensity: 40%
23	28min	Phase 1	5Hz (5mn)	300	EMS	Preparation phase
		Phase 2	55Hz (10s)	300		Intensity: 50% Up: 1.5s Down: 0.75s
		Phase 3	6Hz (8s)	300		Intensity: 25%
		Phase 4	3Hz (10mn)	300		Phase of relaxation Intensity: 40%
24	32min	Phase 1	5Hz (5mn)	300	EMS	Preparation phase
		Phase 2	75Hz (6,3s)	300		Intensity: 50% Up: 1.5s Down: 0.75s
		Phase 3	4Hz (8s)	300		Intensity: 25%
		Phase 4	3Hz (10mn)	300		Phase of relaxation Intensity: 40%

25	28min	Phase 1	5Hz (5mn)	370	EMS	Preparation phase
		Phase 2	55Hz (10s)	370		Intensity: 50% Up: 1.5s Down: 0.75s
		Phase 3	6Hz (8s)	370		Intensity: 25%
		Phase 4	3Hz (10mn)	370		Phase of relaxation Intensity: 40%
26	28min	Phase 1	5Hz (5mn)	370	EMS	Preparation phase
		Phase 2	55Hz (10s)	370		Intensity: 50% Up: 1.5s Down: 0.75s
		Phase 3	6Hz (8s)	370		Intensity: 25%
		Phase 4	3Hz (10mn)	370		Phase of relaxation Intensity: 40%
27	30min	Phase 1	50Hz	200	EMS	Up: 2s Con: 8s Down: 2s Rest: 4s

5.2. DOWNLOADING THE MOBILE APPLICATION

Download the application MyTens your phone or tablet:
- through the Apple store or Google Play Store.

5.3. STARTING AND USING MYTENS / MYTENS PRO.

Check the MyTens App User Manual online before to start using the device

6. ACCESSORIES

MyTens and MyTens pro are supplied with the following unique accessories developed for optimal results:

- Micro-USB cable for charging the device.

- BewellConnect® case.
- MyTens: One pair of electrodes.
- MyTens Pro: two pairs of electrodes.

The electrodes eventually wear out and need to be replaced. It is recommended to replace the electrodes after about 15 to 20 uses. For purchasing information, check BewellConnect® online store or your reseller.

7. STORAGE AND MAINTENANCE

Replace the electrodes in their protective bag and close it well to avoid dust on the electrodes. We recommend storing the electrodes in a place where the temperature is between 41 - 80.6°F / 5 - 27°C. Avoid extreme heat and exposure to direct sunlight. Do not keep the electrodes in the freezer or the refrigerator.

7.1 CLEANING

Warning: Never immerse the product in water and do not rinse it. Never immerse the electrodes in water. Do not use cleaning agents other than those mentioned below. They could cause serious damage to the equipment.

1. Using a damp cloth and mild detergent (eg: dishwashing liquid), clean the device. It is also possible to use isopropyl alcohol at 70° (IPA).
2. If dirt gets on the electrodes, put a drop of water on your finger and gently remove the dust from the surface. Do not use soap or alcohol to clean the electrodes.

8. TROUBLESHOOTING GUIDE

PROBLEM	POSSIBLE CAUSE	SOLUTION
The battery indicator on the device is displayed in orange continuously	Low or discharged battery	Recharge the device
Orange battery light and then the device switches off	Low or discharged battery	Recharge the device

Treatment which seems to be different from or less pleasant than the previous one	1. Electrodes not placed in the correct location 2. Current Intensity too low or very high 3. Expiry date on the electrodes exceeded 4. Worn electrodes 5. Electrodes too dirty 6. You did not correctly apply the electrodes to the skin	1. Switch off the device, take off the electrodes and move them slightly on the skin surface. 2. In the App, modify the intensity level. 3. Remove the electrodes and replace them with new ones. 4. Replace the electrodes. 5. Clean the electrodes (paragraph 5.2 "Cleaning"). If the problem occurs again, replace the electrodes. 6. Switch off the unit, take off the electrodes and make sure to place them properly on the skin.
Unpleasant sensation when touching the connectors of the device or the electrodes	The session has not been stopped in the application or the on/off button has not been pressed before removing the device	Always end the session in the application or remember to press the on/off button before removing the device.
The APP can't connect to the BW-TSX device	The Bluetooth connection failure.	Restart the BW-TSX and the APP.
The APP can't control the BW-TSX device, but the APP can connect the BW-TSX device.	The Bluetooth communication failure.	Restart the APP, and then connect the device again.

More Troubleshooting guide online (Check the FAQ's)

9. TECHNICAL CHARACTERISTICS

UNITS (CONNECTORS)	MyTens: 2 devices (1.9 in x 2.1 in / 49.9 mm x 55.8 mm) MyTens Pro: 4 devices (1.9 in x 2.1 in / 49.9 mm x 55.8 mm)
PROGRAMS	19 programs (11 TENS + 8 EMS programs)
CLASSIFICATION	Medical device Class II
BATTERY	3.7V/500mAh Li-Ion
OUTPUT VOLTAGE	max 60mA ±10% (1000Ω load)

CONDITIONS OF USE	41°F to 104°F / 5°C to 40°C with a relative humidity of 15% to 93%, atmospheric pressure of 700 to 1060 hPa
STORAGE CONDITIONS	14°F to 122°F / -10°C to 50°C with a relative humidity of 10% to 90%, atmospheric pressure of 700 to 1060 hPa
DIMENSIONS	17.3 x 2.1 x 0.5 in / 440 x 54 x 13mm cable included
WEIGHT	0.1lb / 55g (±20%) by unit
IP CLASSIFICATION SYSTEM	IP22
WAVE FORM	Biphasic square wave
IMPULSE WIDTH	150µs~370µs (±10%)
HEART RHYTHM	2Hz~100Hz (±10%)
TREATMENT TIME	1-90 minutes (±10%) given by the application
ADAPTER	1pc (Input : 100-240 V AC 50/60Hz Output: 5 V DC 300mA)

- Service life of electrode pads: 15 to 20 times.
- Applied part: Electrode.
- Size of electrode: 0.98 x 1.10 in / 25 x 28 mm (ellipse).
- Maximum separation distance: 393.7 in / 10 m (In the opened environment).

10. DESCRIPTION OF THE WIRELESS FUNCTION

RF Frequency	2.4GHZ-2.48GHZ
Channel Bandwidth	2 MHz (-20db)
Operation Voltage	DC3.7V
Types of spread spectrum	FHSS
Modulation type	GFSK
Number of channel	40
Chanel spacing	2 MHz
Channel frequency	0-39 Channel 2.402-2.480GHZ
Dwell time (if FHSS)	400ms
Hopping rate (if FHSS)	1600HZ
Antenna gain	0dB

Bluetooth version	V1.0
Maximum Output Powers	0dBm

1. BLE4.0 operate in the unlicensed 2.4 GHz ISM (Industrial Scientific Medical) band. The BLE4.0 wireless technology has 40 channels (3 broadcasting channel, 37 data channel), each channel spacing 2MHz, the entire operating frequency range: 2.402-2.480GHZ. A frequency hopping transceiver is used to combat interference and fading. The way of modulate is GFSK, the output power is 0dbm, and the maximum transmission range is 10m meters.

2. The BW-TSX / TS1 devices use the one-on-one connection of the BLE4.0, and our devices as a peripheral, response to the APP's connection request passively, establish a one-on-one connection with the APP. After the BW-TSX / TS1 devices connect to the APP, the devices will not be connected to other wireless device or controlled. After the APP connect to the peripheral device, the APP will send the command, and when the BW-TSX / TS1 devices ensure this command is right, it will open the control command, then you can control the BW-TSX / TS1 device. Otherwise, the APP can't control the BW-TSX / TS1 device. And make sure that the devices is connected and controlled by our APP.

3. For the reliability of data transmission, the Bluetooth Low Energy 4.0 has take some measures in the radio frequency, baseband protocol, link management protocol to ensure the reliability of the wireless data transmission. And ensure the security of the data with the AES-128 CCM algorithm.

4. This device emits very low levels in the radio frequency (RF) interval and is therefore not liable to cause any interference with nearby electronic equipment (radios, computers, telephones etc.)

5. The BW-TSX / TS1 devices are designed to withstand foreseeable disturbances originating from electrostatic discharges; mains supply magnetic fields or radio frequency transmitters.

6. Despite this, it is not possible to guarantee that the stimulator will not be affected by strong RF (radio frequency) fields emitting from other sources.

7. The BW-TSX / TS1 devices are designed to be used in typical domestic environments approved in accordance with EMC safety standard IEC60601-1-2.

11. IMPORTANT INFORMATION REGARDING ELECTRO MAGNETIC COMPATIBILITY (EMC)

GUIDANCE AND MANUFACTURE'S DECLARATION – ELECTROMAGNETIC EMISSION		
The BW-TSX / TS1 is intended for use in the electromagnetic environment specified below. The customer of the user of the BW-TSX / TS1 should assure that it is used in such an environment.		
EMISSION TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT – GUIDANCE

RF emissions CISPR 11	Group 1	The BW-TSX / TS1 use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The BW-TSX / TS1 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

GUIDANCE AND MANUFACTURE'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The BW-TSX / TS1 is intended for use in the electromagnetic environment specified below. The customer or the user of BW-TSX / TS1 should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.

Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U _T (>95% dip in U _T) for 0.5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 sec	<5% U _T (>95% dip in U _T) for 0.5 cycle 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95% dip in U _T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BW-TSX / TS1 requires continued operation during power mains interruptions, it is recommended that the BW-TSX / TS1 be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

GUIDANCE AND MANUFACTURE'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The BW-TSX / TS1 is intended for use in the electromagnetic environment specified below. The customer or the user of the BW-TSX / TS1 should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
----------------------	-----------------------------	-------------------------	---

Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.6 GHz	3 V RMS outside the ISM band, 6 V RMS in the ISM and amateur radio bands 10 V/m	Portable and mobile RF com- munications equipment should be used no closer to any part of the BW-TSX / TS1 including cables, than the recommended separation distance calculated from the equa- tion applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \text{ 80 MHz to 800 Mhz}$ $d = 2.3 \sqrt{P} \text{ 800 MHz to 2.5 Ghz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
---	---	--	---

			Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol: 
--	--	--	---

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BW-TSX / TS1 is used exceeds the applicable RF compliance level above, the BW-TSX / TS1 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BW-TSX / TS1.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE BW-TSX / TS1

The BW-TSX / TS1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BW-TSX / TS1 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BW-TSX / TS1 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
0.01	150 KHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.1	0.12	0.12	0.23
1	0.38	0.38	0.73
10	1.2	1.2	2.3
100	3.8	3.8	7.3
	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult customer service or an experienced radio/TV technician for help.

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SN:

Retailer's seal Stamp / Sello del revendedor / Cachet du revendeur :

EN: The Bewellconnect® device limited warranty provides that the Bewellconnect® devices will be free from material defects impacting performance in accordance with the documentation for a period of one (1) year when used under normal conditions consistent with the documentation, subject to the terms and limitations set forth below.

The Bewellconnect® device limited warranty is granted only to the original purchaser of the Bewellconnect® device, and is non-transferable. For the sake of clarity, in the case of a refurbishment of the Bewellconnect® device for the use of a different end-user during the original one-year term of the Bewellconnect® device limited warranty, the Bewellconnect® device limited warranty shall continue to be in effect for the remainder of such original one-year term. A purchase invoice or other proof of purchase will be required for after-sales service, in accordance with the Bewellconnect® device limited warranty. The purchase invoice is mandatory to be able to benefit from the Bewellconnect® device limited warranty.

The Bewellconnect® device limited warranty will be invalidated where serial numbers on products are changed, replaced, illegible, absent or if a repair has been made without any outcome by any unauthorized service, including the user.

The Bewellconnect® device limited warranty covers only defective equipment or parts, and only defects that arise during normal use of the product for its intended purpose in accordance with the documentation. It does not cover damage caused during the dispatch or carriage of the device,

caused by repairs made by a distributor, by alterations made, by the connection of equipment not approved by Bewellconnect®, or caused by use that is contrary to the documentation. Furthermore, the Bewellconnect® device limited warranty does not cover damage related to items being dropped, improper handling, improper installation, fire damage, flood, lightning, or any other natural disaster. The Bewellconnect® device limited warranty does not cover product packaging, accessories, or any defective appearance due to the display of the product for sale, in showrooms, retail outlets, for demonstration purposes, etc. Normal use, cleaning and replacement of parts with normal wear and tear are not covered by the terms of this guarantee.

Repair or replacement of the defective Bewellconnect® device is the customer's sole remedy for the Bewellconnect® device limited warranty, or, if the defective Bewellconnect® device cannot be repaired or replaced, refund of the amounts paid by the customer for the defective Bewellconnect® device will be the customer's sole remedy.

BEWELLCONNECT® PROVIDES ONLY THIS BEWELLCONNECT® DEVICE LIMITED WARRANTY, AND MAKES NO OTHER WARRANTIES OF ANY KIND. WITHOUT LIMITING THE FOREGOING, BEWELLCONNECT® EXPRESSLY DISCLAIMS ALL OTHER WARRANTIES, WHETHER ORAL OR WRITTEN, EXPRESSLY IMPLIED, STATUTORY OR OTHERWISE, INCLUDING, WITHOUT LIMITATION, THE IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, AND ALL WARRANTIES ARISING FROM COURSE OF PERFORMANCE, COURSE OF DEALING OR USAGE IN TRADE. BEWELLCONNECT® SHALL NOT BE RESPONSIBLE FOR ANY WARRANTIES OR GUARANTEES GRANTED BY ANY DISTRIBUTOR OR OTHER THIRD PARTY.

TO THE MAXIMUM EXTENT PERMITTED BY APPLICABLE LAW, IN NO EVENT SHALL BEWELLCONNECT® OR ANY OF ITS AFFILIATES BE LIABLE WITH RESPECT TO ANY BEWELLCONNECT® DEVICE FOR (a) ANY SPECIAL, INDIRECT, INCIDENTAL, CONSEQUENTIAL, RELIANCE OR PUNITIVE DAMAGES, INCLUDING, WITHOUT LIMITATION, LOSS OF PROFITS, OR (b) ANY DIRECT DAMAGES IN EXCESS OF THE PRICE THAT CUSTOMER PAID FOR THE BEWELLCONNECT® DEVICE

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EN: The device, accessories, and the packaging must be disposed of correctly at the end of usage. Please follow local ordinances and regulations for disposal.

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EN: *BewellConnect® is a brand registered by Visiomed Group SA. The constructor keeps the right to modify the photographs and characteristics without previous notice. Registered design. The logos and trademarks are the property of their respective owners.*

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