easyScan™

EC REP

VISIONED GROUP SA
112, Avenue Kléber
75116 PARIS - France
contact@visiomed-lab.fr
www.visiomed-group.com
SAV : 0 800 680 080

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Chamberlain, J.M., et al., Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer, *Annals of Emergency Medicine*, January 1995.

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1. Chamberlain, J.M., Terndrup, T.E., New Light on Ear Thermometer Readings, *Contemporary Pediatrics*, March 1994.

2. Chamberlain, J.M., et al., Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer, *Annals of Emergency Medicine*, January 1995.

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By Visiomed

Dear Customer,

You have just purchased EasyScan Pocket infrared thermometer VM-ZXP1, the pocket frontal thermometer by Visiomed. In order to use it optimally and efficiently, we recommend you to carefully read these operating instructions.

The manufacturer reserves the right to modify the product technical specifications without prior notice.

The manufacturer cannot be held responsible for any injury or damages caused as a direct or indirect result of the following cases:

- Improper use of this device,
- Use of inappropriate external components.
- Failure to follow the instructions or misinterpretation of the information presented in this manual.
- Modification or reconfiguration of the device and/or its accessories, without the written consent of the Manufacturer.

TABLE OF CONTENTS

1. SAFETY INSTRUCTIONS.....	21
1.1. WARNINGS.....	21
1.2. CAUTIONS.....	21
1.3. GENERAL PRECAUTIONS	21
2. PRODUCT INTENDED USE	23
3. INTRODUCTION.....	23
4. PRECAUTIONS BEFORE USE.....	24
5. SUGGESTIONS BEFORE MEASURING.....	25
6. MENU SETTING & FUNCTIONS	25
6.1. MAIN FUNCTIONS	25
6.2. CHOICE OF TEMPERATURE UNIT	25
6.3. HOW TO MEASURE TEMPERATURE	25
6.3.1. FOR A FOREHEAD MEASUREMENT.....	25
6.3.2. DIAGNOSTIC OF THE MEASUREMENT.....	26
6.4. DATA MEMORY	27
7. BATTERY REPLACEMENT	27
8. UPKEEP, STORAGE AND CALIBRATION.....	27
9. PRODUCT SPECIFICATION	28
10. TROUBLESHOOTING.....	30

1. SAFETY INSTRUCTIONS

For your safety, to prevent injuries and avoid situations that could damage the device, read carefully all the warnings and precautions in this manual.

Consult your doctor, if you have a serious illness or if you want to know more about using this medical device.

1.1. WARNINGS

- Do not use the unit near the air outlet of air fan or air condition.
- The measurement results are for informational purposes only. If unsure about your results, please contact your physician.
- Do not use the unit if the probe is damaged.
- This product should not encourage you to self-medicate or adapt your treatment.

- Do not touch or blow on the infrared sensor. A dirty or damaged infrared sensor may cause inaccurate results.
- This device have to be used by an adult.

1.2. CAUTIONS

Constraints before taking the temperature:

- Using ice to cool the temperature may cause a lower measuring temperature than its real value.
- If you have high fever or a long time fever, contact your doctor or go to the hospital especially for infants or children.

1.3. GENERAL PRECAUTIONS

- Use this unit only for its intended purpose as described in this manual.
- Follow the maintenance instructions specified in this manual.
- This unit can be used for personal use at home.
- This device is designed for external use.
- Use and storage: Refer to the technical features (range of room temperature, relative humidity, atmospheric pressure).
- This unit should always be stored in a clean, dry place.
- Do not expose this unit to sunlight or water.
- Do not use this unit outdoors.
- Do not expose this unit to electrical shocks.
- Never drop this unit.
- Do not touch the screen with the by finger.
- Do not attempt to open the device. In case of problems, contact the dealer.
- Keep the device out of the reach of children. This contains small pieces that can be swallowed. If swallowed, call a doctor immediately.
- In case of abnormal operation of the device, stop using and contact the After-Sales Service. See the "Troubleshooting Guide" section.
- This electrical medical equipment requires special precautions regarding electromagnetic compatibility. It must be installed and used in accordance with electromagnetic information provided in the EMC Data table.
- Readings may be distorted if the unit is used near a television, microwave oven or any other device with an electrical field.

2. PRODUCT INTENDED USE

EasyScan Pocket, connected Infrared Thermometer is intended for human body temperature measuring for all ages (for infant, children and adult). Body temperature can be taken from forehead. It is intended to be used for family healthcare and hospital clinical reference.

3. INTRODUCTION

The normal body temperature is located within a range, and a fixed value. Everyone has their own temperature range. When the body temperature rises or falls, the body initially regulates its temperature from the brain. Since the forehead and the ear are the body parts nearest to the brain, they are the first to be informed of a temperature rise.

Practical considerations when taking a temperature:

- In the interests of precision and measurement quality, it is essential that each user is adequately informed and trained in the temperature measurement technique when using this device.
- It is essential to remember that the procedure of taking a temperature, although simple, should not be ignored.
- The patient should not carry out strenuous activities before taking their temperature.
- Temperature fluctuations due to physiological factors must be taken into consideration in the evaluation of the results:
 - The temperature increases by 0.5°C/0.9°F between 6:00 and 15:00,
 - The average temperature of women is about 0.2°C/0.36°F and varies according to their ovarian cycle,
 - It increases by 0.5°C/0.9°F during the first half of the cycle and during the first months of pregnancy.
 - Age, sex, type and thickness of skin.
- The temperature of a seated person is lower by around 0.3°C to 0.4°C/32.54°F to 32.72°F than that of a person standing up.

Haut Conseil de Sante Publique [online]. Health Recommendations of the National Heatwave Plan 2014. P123 [Consulted on 25 July 2016] Available at https://www.hcsp.fr/Explore.cgi/Accueil/Recommendations_sanitaires_Plan_canicule_2014.Pdf

Normal temperature by age

Age	°C	°F
0-2	36.4 - 38.0	97.5 - 100.4
3-10	36.1 - 37.8	97.0 - 100
11-65	35.9 - 37.6	96.6 - 99.7
> 65	35.8 - 37.5	96.4 - 99.5

Chamberlain, J.M., et al., Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer, Annals of Emergency Medicine, January 1995.

Normal temperature according to procedure

Measurement method	Normal T°
Rectal ⁽¹⁾	36,6 – 38°C / 97,88°F – 100,4°F
Oral ⁽¹⁾	35.5 – 37°C / 95,9 – 99,5°F
Axillary ⁽¹⁾	36,5 – 37,5°C / 97,8°F – 99,5°F
Aural ⁽²⁾	35,8 – 38°C / 96,44°F – 100,4°F

1. Chamberlain, J.M., Terndrup, T.E., New Light on Ear Thermometer Readings, Contemporary Pediatrics, March 1994.

2. Chamberlain, J.M., et al., Determination of Normal Ear Temperature with an Infrared Emission Detection Thermometer, Annals of Emergency Medicine, January 1995.

4. PRECAUTIONS BEFORE USE

EasyScan Pocket is pre-configured at the factory. It is not necessary to calibrate the device for the first use.

To obtain reliable, stable results, please follow below instruction.

- Check the unit for any damage in appearance, and check the LCD screen full display as below picture. Do not use the unit if any abnormality and contact your retailer. If battery capacity is low, follow the instruction to change to a new battery.

- Before any new temperature measurement, wait for the screen to turn off.

- It is advisable, with each change of environment with a significant difference in ambient temperature, leave the EasyScan Pocket at this ambient temperature for at least 30 minutes before use. It is important

to wait a period of 1 minute between measurements. If the person come from the different ambient temperature than the unit, please stay at least 5 minutes before taking temperature.

5. SUGGESTIONS BEFORE MEASURING

1. The measurement should be taken at rest, to yield accurate results. Body temperature is likely to increase just after exercise or after a bath.
2. Before taking a temperature, in order to obtain a stable, reliable result:
 - Make sure to brush away the hair.
 - Wipe off sweat.
 - Avoid air flows (e.g. supplemental oxygen, air conditioning, etc.).
 - Wait a period of 1 minute between each measurements.

6. MENU SETTING & FUNCTIONS

6.1. MAIN FUNCTIONS

1. Display of the temperature reading in degrees Celsius or Fahrenheit: Temperature range from 30°C to 43°C (86°F to 109.4°F).
2. Memorization of the last 10 readings taken.
3. Audible alarm during temperature measurement based on the result.

6.2. CHOICE OF TEMPERATURE UNIT

Open the battery cover, press the C°/F button to select the temperature unit. The display shows C° to display degrees Celsius, °F to display degrees of Fahrenheit.

6.3. HOW TO MEASURE TEMPERATURE

Insert the batteries. For the first use or when inserting new batteries, wait about 10 minutes to preheat the device.

6.3.1. FOR A FOREHEAD MEASUREMENT

Picture A

1. Press the measurement button once and make sure the logo ☺ of the front temperature sensor is displayed on the screen.
2. Point the infrared sensor towards the forehead (picture A), it should touch the skin slightly. Make sure

to stay still during the measurement.

3. Press the measuring button to start for one second. The unit will beep once, it indicates the measurement is completed, and the LCD screen will display the measured temperature value (refer to the diagnostic table for measurement below).
4. The result is automatically saved. In the case of consecutive measurements of temperature, the results may not be identical. Wait at least 1 minute between each measurement.
5. If the temperature exceeds the warning temperature value, a beep sounds four times.
6. Holding the device in the hand, during a long time, will make the detector read an ambient temperature higher than the real one. That may lead to an inaccurate measurement.

6.3.2. DIAGNOSTIC OF THE MEASUREMENT

T (°C)	Beep	Diagnostic
30 ≤ T° < 35	Di-Di-Di-Di	Hypothermia ¹
35 ≤ T° ≤ 36,2	Di	Temp. to watch
36,3 ≤ T° ≤ 37,3	Di	Normal temp. ²
37,4 ≤ T° ≤ 37,9	Di	Temp. to watch
38 ≤ T° ≤ 43	Di-Di-Di-Di	Fever ³

1. REULER, James B. Hypothermia: pathophysiology, clinical settings, and management. *Annals of internal medicine*, 1978, vol. 89, no 4, p. 519-527.

2. Koushik Maharatna,Silvio Bonfiglio, Systems Design for Remote Healthcare, Springer Science & Business Media, 13 nov. 2013 - 332 pages.

3. Wunderlich, C.A. and J.C. Reeve, The course of the temperature in diseases: a guide to clinical thermometry. *American Journal of Medical Science*, 1869, 57: p. 423-47 . Various sources have slightly different temperatures for fever.

When the temperature is <30 ° C, Lo is displayed on the screen.

When the temperature is > 43 ° C, Hi is displayed on the screen.

6.4. DATA MEMORY

Open the battery compartment to access the setting buttons. Press the MEM button, the LCD screen reads and displays the last temperature reading.


1. " --- " Will be displayed if no data is stored
2. Press the MEM button to scroll through the previous data, then press the °C/F to scroll through the next data. 10 measurements are memorized.

Note: When changing batteries, the stored data will be erased.

7. BATTERY REPLACEMENT

Remove the battery cover and use a tool to take out the drained battery (below picture C), then place a new battery, taking care to respect the polarity. The LCD screen will display all the symbols on the screen.

Use only new batteries (1 x CR2032), never use rechargeable batteries. Replace the cover after inserting the battery (below picture D).

IMPORTANT:

- Respect the polarity of the batteries. A polarity error can cause damage and void the warranty of your unit.
- Do not use battery other than the battery mentioned, do not recharge non-rechargeable battery, do not incinerate them.
- Remove the battery when not using the thermometer over a prolonged period. This may damage the product.

8. UPKEEP, STORAGE AND CALIBRATION

UPKEEP:

1. Remove any spots from the unit with a soft, dry cloth. For stubborn spots, clean the unit with a wrung cloth after soaking in water or a neutral solvent, then wipe with a soft, dry cloth.
2. Clean the infrared sensor as follows: wipe the surface thoroughly with a

- pad, a soft cloth slightly moistened with 70 ° alcohol. The thermometer can be used once the alcohol is completely evaporated.
3. Do not clean the infrared sensor with a tissue or paper towel that may cause scratches and therefore incorrect results.
 4. Do not use other chemical reagents to clean the lens.
 5. Do not immerse the device in any liquid.
 6. Do not clean the display of the device with alcohol.
 7. Do not use benzene, thinner, gasoline or other strong solvents to clean the unit.
 8. Do not use ultrasonic cleaning to clean the unit.
 9. To disinfect, use a cotton swab or a soft dry cloth soaked in 70% alcohol to clean the unit.
 10. If you do not use the device for a long time, it is recommended to clean the unit one time each month.

STORAGE:

1. Do not expose the thermometer to direct sunlight, extreme temperatures conditions ($T > 60^\circ\text{C}$ or $T < -20^\circ\text{C}$), or under conditions of high humidity (above or equal 95%), dust, near fire or high concentrations of salt in the air; do not subject it to inclinations, vibrations or direct impacts; Preservation of chemicals or the presence of corrosive gases.
2. Remove the batteries from the housing if you do not plan to use the unit over a long period.

CALIBRATION:

This thermometer is pre-configured at the factory. If the operating instructions are followed, a re-calibration is not necessary. If you have doubts about the accuracy of temperature readings, please contact your retailer.

Do not disassemble the thermometer yourself - this will void your warranty.

9. PRODUCT SPECIFICATION

Product name	EasyScan Pocket Infrared Thermometer
---------------------	--------------------------------------

Model	VM-ZXP1
Classification	Class IIa
Product size	76mm × 50,6mm × 25,5mm (L×W×H)
Measuring range and accuracy	30°C ~ 34,9°C : ±0,3°C 35,0°C ~ 42°C ±0,2°C 42,1°C ~ 43°C : ±0,3°C
Repeatability	≤ ±0,2°C/F
IR wavelength	5 to 14μm
Emissivity	0,95
Automatic shutoff	15s
Data storage	Last10 readings
Displaying of reading	°C/F
Display resolution	32.1°F / 0.1°C
Internal time of each measurement	>60s
Battery	CR2032×1
Voltage	2,5V ~ 3,3V
Dynamic current	≤ 15μA
Quiescent current	≤25mA
Switch current	≤ 10mA
Consumption	≤ 90mW
Product weight	38g
Normal condition of use	Temperature: 15°C ~ 40°C Humidity: 10 % ~ 85 %RH (non-condensing) Air pressure: 70kPa ~ 106kPa

Storage conditions	Temperature: -20°C ~ 60°C Humidity: ≤ 95 %RH Air pressure: 50kPa ~ 106kPa
Water proof	IP22
Warning temperature	< 35 °C ≥ 38 °C
Expected life time	5 years

The screen display Er2	Ambiant temperature is higher than 40°C	
------------------------	-----------------------------------------	--

10. TROUBLESHOOTING

Problem	Possible cause	Solution
The device does not switch on	The batteries are discharged	Put new batteries in the device
	The batteries are not properly installed	Install the batteries properly, following the correct polarity
Inconsistent measure	Temperature reading hampered by hair, sweat,...	Make sure there is no obstruction prior to taking the temperature.
	Temperature reading hampered by an air flow...	Be careful not to be in the axis of an airflow, this could cause interference with the infrared system.
The display indicates temperatures above 95°	The temperature is indicated in Fahrenheit.	Change the unit of measure to Celsius.
The screen displa Er1	Ambiant temperature is lower than 15°C	

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1. REULER, James B. Hypothermia: pathophysiology, clinical settings, and management. Annals of internal medicine, 1978, vol. 89, no 4, p. 519-527.

2. Koushik Maharatna,Silvio Bonfiglio, Systems Design for Remote Healthcare, Springer Science & Business

Media, 13 nov. 2013 - 332 pages.

3. Wunderlich, C.A. and J.C. Reeve, The course of the temperature in diseases: a guide to clinical thermometry. American Journal of Medical Science, 1869. 57: p. 423-47. Various sources have slightly different temperatures for fever.

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EMC DATA

EMC data**Guidance and manufacturer's declaration – electromagnetic emission - for all EQUIPMENT AND SYSTEMS**

The model VM-ZXP1 is intended for use in the electromagnetic environments specified below. the customer or the user of the model VM-ZXP1 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV, +/ - 8 kV contact +/- 2 kV, +/ - 4 kV, +/ - 8 kV +/- 15 kV air	+/- 6 kV, +/ - 8 kV contact +/ - 2 kV, +/- 4 kV, +/ - 8 kV +/- 15 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	+/- 2 kV for power supply lines +/- 1 kV for input/output lines	Not applicable	Not applicable
Surge IEC 61000-4-5	+/- 1 kV line to line +/- 2 kV line to earth	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_1 (>95% dip in U_1) for 0.5 cycle <5% U_1 (>95% dip in U_1) for 1 cycle >70% U_1 (30% dip in U_1) for 25/30 cycles <5% U_1 (>95% dip in U_1) for 5/6 sec	Not applicable	Not applicable
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m, 30 A/m	3 A/m, 30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment

NOTE: U_1 is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic emission - for all EQUIPMENT AND SYSTEMS

The model VM-ZXP1 is intended for use in the electromagnetic environments specified below. the customer or the user of the model VM-ZXP1 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
	3 Vrms 150 kHz to 80 MHz		Portable and mobile RF communications equipment should be used no closer to any part of the models model name, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
	6 Vrms in ISM bands	Not applicable	Recommended separation distance
	3V/m, 20 V/m 80 MHz to 2.5 GHz	Not applicable	$d=[3,5/V_1] \times P^{1/2}$
Conducted RF IEC 61000-4-6	10 V/m 80 MHz to 2.7 GHz	3V/m, 20 V/m 80 MHz to 2.5GHz	$d=1.2 \times P^{1/2} 80 \text{ MHz to } 800 \text{ MHz}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7GHz 385MHz-5785MHz	10 V/m 80 MHz to 2.7GHz 385MHz-5785MHz	$d=2.3 \times P^{1/2} 800 \text{ MHz to } 2.7 \text{ GHz}$
		Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment (Refer to table 9 of IEC 60601-1-2:2014)	where P is the maximum output power rating of the transmitter In watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,a should be less than the compliance level in each frequency range.b Interference may occur in the vicinity of equipment marked with the following symbol: 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the models model name are used exceeds the applicable RF compliance level above, the models model name should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the models model name.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

Recommended separation distances between portable and mobile RF communications equipment and the models VM-ZXP1

The models **VM-ZXP1** are intended for use in electromagnetic environment in which radiated RF disturbances is controlled. The customer or the user of the models **VM-ZXP1** can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the models **VM-ZXP1** as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=1.2 \times P^{1/2}$	80 MHz to 800 MHz $d=1.2 \times P^{1/2}$	800 MHz to 2.5GHz $d=2.3 \times P^{1/2}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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EN: Visiomed® will repair or replace this product free of charge in the case of defective parts or manufacturing defects, in accordance with the conditions mentioned below as follows:

DURATION: 24 MONTHS RETURN TO WORKSHOP

LIMITS AND EXCLUSIONS: This guarantee concerns only the original final purchaser. A purchase invoice, or another proof of purchase, with this guarantee card will be required to obtain an after-sales service, in accordance with this guarantee. This guarantee card will not be extended to another person only the original final purchaser. This guarantee becomes void if the serial numbers on the product are modified, replaced, illegible, absent, or if repair has been carried out by a service not approved, including the user.

This guarantee covers only the defects of the material or parts, occurring during normal use of the product. It does not cover the damage caused

during the transport of the apparatus, causes due to repairs being carried out by the distributor, by any modifications undertaken, any connection of equipment not approved by Visiomed®, or causes contrary to those written in the user manual or notice. Moreover, the present guarantee does not cover damage due to falls, bad handling, bad installations, damage by fire, floods, lightning, or any other natural disaster. This guarantee does not cover the packing of the material, the accessories, the defects caused by commercial exposure of the product, show room, sale space, demonstration etc... Normal maintenance, cleaning and the replacement of parts where wear is normal, are not covered by the terms of this guarantee. Visiomed® and its representatives and agents will not in any case be held responsible for any damage and consecutive damages due to the mishandling of this product. This guarantee is the only valid one at Visiomed®, any other guarantee (commercial guarantee) except this one will not be taken into account.

IMPORTANT: During the guarantee period if you are dissatisfied with the repairs of this product, please contact the Visiomed® customer service.

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EN: The symbol on the product or on its packaging indicates that this product must not be disposed of with your other household waste. Instead, it is your responsibility to dispose of your waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local city office or your household waste disposal service.

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